COMPACT-TENS ECO-DE®

User Manual

Index

Mechanism	2
Features	2
Instructions for correct operation.....	2
Usage	3
Safety Precautions.....	4
Maintenance and storage	4
Trouble Shooting	5
Charging tips.....	5

Mechanism

Based on physics, bionics, bio-electricity, human engineering principles and combined with modern microelectronic technology, research and development of high-tech, Compact-Tens is a sport body sculptor device. It is controlled by the microcomputer processor, equipped with six groups of different methods of automatic procedures for body sculpting. It can produce different frequencies and intensities of bio-pulse; these pulses will stimulate the safety of the body's motor neurons, so that muscle contraction exercise and burn fat in the layer of oxygen-rich environment. Bio-pulse signals effectively breaks down fats, consumption of excess calories.

Many professional athletes are using this advance technology to stay in shape. Now, you can sit comfortably at home to enjoy the benefits of body sculpting method. Put the pads on the parts of your body you want to exercise, and after a period, you will have a perfect shape. Moreover, it can eliminate back pain, neck stiffness and other diseases of modern civilization.

Features

1. Without side effects, it does not produce dependence.
2. Modern design and 6 different sets of body sculpting modes (Intelligent mode).
3. Gentle-start design
4. 10 strength can be adjusted
5. Built-in timer, sports instrument off automatically after 15 minutes of work.
6. Standby mode, 20 seconds after motion instrument (it will automatically shut down)

Instructions for correct operation

This button is used to select the mode; six kinds of preset mode, the corresponding button will be lighted.

ON/: turn on, increase output strength. Press this button once at the beginning, and then press again to increase intensity.

OFF/: press this button once, corresponding to decrease output strength, until turning off the device.

LED I: Sport fat burn patterns

LED II: Deep fat burn patterns

LED III: Shock shape pattern

LED IV: Massage mode

Usage

1. Before using, consume the rest of electricity. Then charge.
2. Aim magnetic clasp of the patch at the magnetic buckle on the bottom of the main machine.
3. Tear off the smooth protective film of the patch. Do not throw away, keep it.
4. Put the device together with the patch on the part of your body you want to exercise. Make sure it makes contact to your body.
5. Press ON
6. After the light, press MODE to select function.
7. After selecting the mode, press ON button to increase intensity.
8. To reduce intensity press OFF button.
9. To switch the device off, press OFF button until it stops. It has an automatic shut down after 15 minutes.
10. Put the film back into the patch after using and storage carefully.

Safety Precautions

Never use the equipment with the following instruments:

- 1) Heart rate regulator and body built-in electronic therapeutic instruments.
- 2) Artificial cardio pneumatic pacemaker, electronic life support and treatment instruments.
- 3) Cardio meter and wearing-on type electronic therapeutic instruments.

Malfunction or danger may occur for the patients who have been using the above listed instruments together with our equipment.

Follow doctor's advice for the patients who are under medical treatment or for sufferers from the following:

- 1) Malignant tumour
- 2) Menstrual phase or pregnancy
- 3) Do not adopt the equipment near of heart disease patients
- 4) Fever above 38°
- 5) Abnormal cutaneous anaphylaxis
- 6) Do not adopt the equipment in the eye, inside the mouth, genital parts and external skin trauma.
- 7) Do not use the equipment simultaneously with other instruments or coated medicine.

Maintenance and storage

- 1) Clean the main part of the device with a soft cloth. Never use thinner gasoline or other toxic liquids.
- 2) Cleaning the electrode: a contaminated electrode may be cleaned with wet cloth or rinsed with fresh water, it can be used again after dry.
- 3) Storage: place the device out of children's reach; keep it away from moisture, high temperatures, direct sunlight or water, etc. Take batteries out of the device when it is not in use for a long time.

Trouble Shooting

Status of the Product	Cause	Solutions
Cannot turn on	Battery used up	Check the cause and charge the battery
	Adhesive sides are damaged	Replace with new pads
Stimulation becomes weak or tingling	Electrode contacting skin correctly	Stick the mat firmly onto the skin
	Dirt on mat surface	Clean or change mat
	Intensity too weak	Adjust ON button or change the battery
Treated surface flushing, itching or tingling	Skin reaction against mat contact	Check skin allergies. For slight or severe allergy, stop treatment
	Therapy is too long	Shorten time
	The adhesive is dirty or dry	Wash it with water, use it after dry

Charging tips

1. When charging, the red light will be on.
2. After red light turns green, means that the charging is complete, unplug the charger and charging connector.
3. When the green light is on, unplug the charging plug, to avoid prolonged transition charge.

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WARRANTY

CONDITIONS:

- A. This product guarantee will be accepted by ECOGEST TRADE S.L., provided that it is correctly completed and the attached card has been stamped and send to ECOGEST TRADE S.L.
- B. The responsibility of ECOGEST TRADE S.L. C/ En Proyecto 14, 46190 Ribarroja de Turia, Valencia (Spain) it is limited to the repair of the articles that has some deficiency of technical character or of manufacture. The return postage will be borne by the seller.
- C. All repairs which are carried out under the present guarantee should be performed by an establishment which is duly authorized by ECOGEST TRADE S.L.
- D. This guarantee will not be considered valid in those cases which result from the misuse.
- E. Incorrect application or manipulation of the product.
- F. Validity of guarantee for 2 years from the delivery date.

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