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Energía: <3W

ELECTROSTIMULATOR MASSAGER **ECO-308B**

Multi-function electrostimulator massager with double channel.

Through electrical stimulation is accomplished separating the fats that are deposited in certain areas to exert a kind of drainage in order to eliminate and prevent them from falling back into the body.

By removing the fat, helps muscles to assert, reactivate the bloodstream to be carried through the blood toxins.

Sessions last between 30 and 45 minutes depending on the area treated and the results you seek, for example, if you only want to reduce cellulite and toning, buttocks, time varies, and the frequency of sessions.

Electrostimulation benefits:

- Acts directly on problem areas where fat deposits are located and therefore the results are more evident.
- You can lose inches and kilos.
- Help as a treatment for cellulite.
- Defines and helps improve muscle power.
- In addition to undo the fat deposits, helps expel body and firm muscles.
- There is a risky treatment does not require anesthesia, painless and safe
- It helps to regenerate and oxygenate their tissues

Contraindications of electrostimulation:

- This method is not recommended for people with epilepsy, pacemakers, tumors, varicose veins, thrombosis, bleeding, infections, pregnant women with fever, emotionally disturbed or under 10 years, etc

Parts names:

- | | |
|---------------------------------|---------------------------|
| 1. ON/OFF | 12. WAIST |
| 2. INTENSITY SWITCH (CHANNEL 1) | 13. LEG |
| 3. ARM | 14. TIME + |
| 4. SLIMMING | 15. TIME - |
| 5. FITNESS | 16. SHOULDER |
| 6. FEET | 17. INTENSITY CHANNEL 2 |
| 7. JOINT | 18. MESSAGE SLIPPER |
| 8. PRESS | (AVAILABLE) |
| 9. KNEAD | 19. ELECTRODE PADS |
| 10. VIBRATE | 20. 2 IN 1 ELECTRODE WIRE |
| 11. THUMP | 21. POWER |

Operating instructions:

- . Open the box and put 4 pieces of 1.5v batteries.
Close. 6 V adaptator can also be applied
- . Stick the electrode pads onto appropriate parts of the body and then connect conducting wire with the electrodes correctly.
- . Open the on-off button of the main machine, and select of the working mode. You can choose different massage modes (press, knead, vibrate and thump) according to your personal preference.
- . The massage speed and intensity can be adjusted flexibly. You can tune the intensity knob of the main machine from lower to higher.
- . You can also press the button "arm", "joint", "feet", "shoulder", "waist" and "leg" to choose different parts of body massage.
- . Default working time of the machine is 30 minutes. The machine also has 5 minutes of self checking control. The machine will cut off once removing the electrode pads from the skin within five seconds automatically
- . It's double channel system can make 2 persons use it at the same time. Intensity and time can be adjusted separately.

. When you stop the machine, please turn the power off

Notes:

. Never use the machine with the following medical appliances:

. pulse regulator

. artificial cardiac pacemaker

. electrocardiograph machine

. Don't put the electrodes onto the point or precordial region, carotid and the waist of pregnant woman.

. Seriously sick person should consult doctors before using the machine

. Don't touch the electrode pads with metal

. Don't use the machine when you feel uncomfortable or if you have skin irritation

. Don't use the machine under wet conditions, such as bathrooms

. never throw the batteries into the fire.

Storage / maintenance:

- . Use wet or neutral cleaner to wipe the machine.
 - . Don't use volatile liquid to clean the machine, such as benzene, thinner or gasoline.
 - . Keep away from magnetic field
 - . Keep the machine away from fire, heat and water
 - . Don't leave the machine in direct sunlight for a long period time
- , Please take out the batteries if the machine is not in use for a long time.

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WARRANTY

CONDITIONS:

- A. This product guarantee will be accepted by ECOGEST TRADE S.L., provided that it is correctly completed and the attached card has been stamped and send to ECOGEST TRADE S.L.
- B. The responsibility of ECOGEST TRADE S.L. C/ En Proyecto 14, 46190 Ribarroja de Turia, Valencia (Spain) it is limited to the repair of the articles that has some deficiency of technical character or of manufacture. The return postage will be borne by the seller.
- C. All repairs which are carried out under the present guarantee should be performed by an establishment which is duly authorized by ECOGEST TRADE S.L.
- D. This guarantee will not be considered valid in those cases which result from the misuse.
- E. Incorrect application or manipulation of the product.
- F. Validity of guarantee for 2 years from the delivery date.

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